



## DEKO 25

The DEKO 25 has been designed to bring reliability, outstanding cleaning and disinfection performance and quality into one compact and economic washer-disinfector.

The DEKO 25 can be situated on a dedicated DEKO stand accessory or a cabinet, as well as, be positioned undercounter.

The DEKO 25 meets the performance and design requirements stated in EN ISO 15883 parts 1 & 2, as well as the safety requirements stated in IEC 61010-2-040.

Utilising the power of the 260 litres per minute circulation pump, process and cleaning results are guaranteed throughout the wash chamber. The capacity also allows for the insertion of additional wash racks to facilitate the washing of cannulated instruments, dental handpieces, hoses, glassware and rigid scopes.

The DEKO 25 comes as standard with pre-set programs installed and protected, however additional tailored programs can be added to meet any process requirements such as dental or laboratory use.

The washing spray patterns are delivered through two rotating spray arms or utilising jet rack connectors and guarantee water contact with all surfaces and crevices of instruments with the most intricate design.



A base cabinet or a stand are available as an accessory.

## LOGIC CONTROLLED

The DEKO 25 presents the latest technology in user interface applications in washer-disinfectors. Information & control panel provides the means to view cycle information real time and view the remaining cycle time. Also included is an automatic cycle counter and the intelligent software provides advice on service actions.

## EASE OF USE

The DEKO 25's door opens & closes smoothly and is locked securely throughout the cycle by the automatic locking system. The clear text on the display and the cycle complete indicator will inform the operator when the cycle has been completed.

The in-built condenser prevents the escape of excessive steam and minimises the moisture within the chamber at the end of the cycle.



## TRACEABILITY AND VERIFICATION

In order to guarantee that the values of the critical process variables were truly obtained and as programmed, the DEKO 25 is equipped with an integrated process verification printer.

The critical process values i.e. temperatures, times, pressures, detergent dosage time and volume, faults, as well as descriptive cycle data can be independently recorded by for your Quality Assurance System.

Archive information includes: institution name, date and time and machine identification, cycle number, selected cycle name, phase descriptions, detergent dose volume and stage, chamber and wash temperature throughout the process and faults.

## DURABLE QUALITY

The DEKO 25 is manufactured using quality polished stainless AISI 304 steel, it is easy to keep clean and is designed to withstand continuous use. Galvanised and painted metals and other non-compatible materials are totally avoided in the manufacture of the DEKO 25. Special attention has been given to sound and heat emissions by using the most efficient insulation and environmentally friendly materials.

Durability and reliability will allow continuous use in any busy and demanding environment, these types of heavy workload have been the key factors in the design of the DEKO 25.





## RELIABLE

The DEKO 25 occupies its space efficiently - one unit requires only W600xD500 mm of floor space. All components are accessible by the quick removal of dedicated service panels and can be made without moving the washer. Efficient diagnostic is achieved by interrogation of the PLC through the simplistic service mode.

Water is introduced directly into the chamber through a type "A" air gap to avoid the need of cleaning and disinfecting separate water holding tanks.

Any information, personnel training and service calls are attended to with the utmost dedication and urgency.

## SECURE

The advanced electronic control system can provide practically an unlimited number of tailored cycles with the 3 most used cycles available at a touch of a button and additional ones by scrolling. Process parameters are protected by security codes. The machine is facilitated by a 3 level access code to prevent unauthorised use of the machine.

Protection to thermo-labile medical devices is provided by a double overheat protection.

To achieve the outstanding, consistent cleaning efficiency the washer disinfector monitors the detergent volume very closely as well as preventing the washer to start if there is insufficient detergent volume to complete the cycle.



## SUMMARY OF FEATURES

- overall dim. W460xD580xH845 mm
- usable wash chamber dim. W265xD390xH390 mm
- powerful circulation water pump, output 260 L/min
- framework, panels and wash chamber all made of stainless steel
- interlocking door
- fully automatic PLC control including back-up batteries for full memory protection
- 3+1 cycle selection buttons (1 scrolling)
- phase indication as text
- remaining cycle time indication
- cycle complete indication as text, a green led and an audible indicator
- fault indication
- multi-level access code system. Special codes for: user, maintenance, authorized person
- practically unlimited number of pre-set cycle library
- customized cycle files available
- access to regulation of critical process parameter values with authorized person entry code only
- cycle in progress: yellow led
- cycle complete: green led
- fault indication: red led
- special service and maintenance control system
  - malfunction identification
  - adjustment of controls
  - calibration of measuring instruments
  - manual override step by step
  - component test mode
  - fault listing
  - updates by a PC
- 1 sliding wash wire basket
- automatically connecting / closing wash arm and rack connectors
- 2 rotating wash arms with removable caps facilitating easy cleaning of the inner surfaces of the wash arm. No tools required for removing the wash arms
- detergent flow control
- adjustable automatic detergent dosing
- indication for insufficient volume of detergent to complete cycle
- boss entry port for test temperature sensor leads
- load temperature protection and overtemperature cut-outs
- automatic foam and washing water pressure control system
- opaque cold and hot water connection hoses 1/2" BSP
- drain outlet hose diameter 20mm
- sound and heat insulated
- height adjustable feet
- back syphon Type A air gap protection

## OPTIONS

- upper basket incl. wash arm
- max 1 extra detergent pump
- max 1 extra adjustable rinse aid pump
- solenoid valve for connection to purified water supply
- cabinet or a stand for process chemical containers (accessory)


## TECHNICAL DATA

Model	DEKO 25
<b>Outer dimensions mm:</b>	
Width	460
Depth	580
Height	845
Weight kg	60
<b>Wash chamber volume L</b>	
	40
<b>Usable inner dimensions mm:</b>	
Chamber height	390
Chamber width	265
Chamber depth	390
Wash levels	1(+1)
<b>Usable height on lower level, mm</b>	
	250
<b>Usable height on upper level, mm</b>	
	118/138
<b>Wash arms max. (2 as std.)</b>	
	3
<b>Auto-connection for injection tube racks</b>	
	1
<b>Loading basket 370x370x80 mm</b>	
	1(+1)
<b>Water volume / one cycle phase, L</b>	
	6.5
<b>Electricity</b>	
240V 1N 50Hz 3.5kW	16A
230V 1N 50Hz 3kW (United Kingdom)*	13A
200V 1N 50/60Hz 3kW*	16A
Please specify and contact for further details.	
<b>Capacity &amp; effects:</b>	
Heating power kW	3
Heating power kW*	2.5
Circulation pump 260L/min, W	540
Outlet pump 28L/min, W	85
Outlet size diameter in mm	20
Outlet hose length 2 m	
<b>Cold water hose L=1,5 m</b>	
	R1/2"
<b>Hot water hose L=1,5 m</b>	
	R1/2"
<b>Purified water connection</b>	
	R3/4"
<b>Min/Max water pressure bar</b>	
	1/8

Subject to changes without further notice. Please ask the supplier for confirmation as required.










<b>Included as standard:</b>	<b>Tests / Approvals:</b>	SIS TR-3:2001
• 1 wash basket	ISO 9001	EN ISO 15883-1, 15883-2
• 2 wash arms	ISO 14001	IEC 61010-2-040
• 1 detergent pump	ISO 13485	<b>CE 0537</b>


Installation example




# ACCESSORIES

1. 381326  
**Top wash basket incl. a wash arm**  
- dim. 380x380x90 mm
2. 374085  
**Lab wash base** - dim. 390x390 mm
3. 381210  
**Instrument tray, 340x250x70 mm**  
381220  
**Instrument tray, 250x170x70 mm**  
381211  
**Tray lid, 340x250 mm**  
381221  
**Tray lid, 250x170 mm**
4. 381215  
**Instrument tray with end plates**  
340x250x70 mm  
381225  
**Instrument tray with end plates**  
250x170x70 mm  
381211  
**Tray lid, 340x250 mm**  
381221  
**Tray lid, 250x170 mm**
5. 381212  
**Fine mesh tray, 330x240x50 mm**  
381222  
**Fine mesh tray, 240x160x50 mm**  
381213  
**Tray lid, 330x240 mm**  
381223  
**Tray lid, 240x160 mm**
6. 6381250  
**Base basket, one included as standard**  
- dim. 380x380x90 mm
7. 7381170  
**Utensil basket with a lid, 340x250x200 mm**
8. 381685  
**Rack for 6 kidney bowls, 330x190x160 mm**
9. 381525  
**Peg rack, 9xH150 mm pegs**  
- dim. 350x210 mm
10. 381214  
**Holder for 2 instrument trays**  
- dim. 320x130x130 mm
11. 381635  
**Tray holder, full size, 360x360x130 mm**  
381640  
**Stainless steel tray, 280x180x16 mm**
12. 381540  
**Container for microscope slides, incl. a lid**  
- dim. 250x170x100 mm






13




19



14




20




15




21



16




22




17



23



18



24

13. 381325  
**Injector rack for dental equipment e.g. angle pieces, turbines for outer cleaning**

14. 381490  
**Injector rack for dental equipment e.g. angle pieces, turbines for outer cleaning**

15. 381491  
**Injector rack for hollow devices eg. needles**

16. 381492  
**Insert for 'Oras' faucet handles and hoses**

17. 381493  
**Rack for rigid scopes**

18. 381494  
**Rack for AN-equipment**

19. 381560  
**Injector rack**  
- 18 jet tubes H150 mm, Ø 8 mm  
- 12 jet tubes H100 mm, Ø 6 mm

381725  
**Jet tube tensioning key**  
N23658 Bottle holder for dia. 8 mm jet tube  
N23659 Bottle holder for dia. 6 mm jet tube

381556  
**Bottle support 4 pegs, for 8 mm jet tube**

381559  
**Bottle support 4 pegs, for 6 mm jet tube**

20. 381218  
**Mirror compartment**  
- dim. 370x185x110 mm

21. 381551  
**Rack for 18 Petri dishes**  
- dimensioned to diameter  
- std dim. 330x300x90 mm

22. 381260  
**Fine mesh tray, 370x370x40 mm**

381536  
**Pressure lid, 370x370 mm**

23. 381475  
**Half-circle basket, D 380 mm, H 135 mm**  
- includes a lid and compartment divider

24. 381227  
**Storage basket for angle pieces and turbines**  
- 30 places, dim. 250x170x70 mm

381230  
**Small box for needles and bore bits**  
- dim. 90x90x40 mm

381235  
**Instrument tray, 250x100x70 mm**

**Industrial, Catering and Medical washing equipment**

72 -74 Woodlands Drive, Braeside, Victoria, 3195  
PO Box 605 Braeside BC 3195

Tel: + 61 3 8586 5444

Fax: + 61 3 8586 5454

info@rhima.com.au

www.rhima.com.au

Rhima Australia Pty Ltd

ABN: 63 007 242 604

